

## A100D(LCD) Portable Fetal Doppler

A100D(LCD) Portable Fetal Doppler is a high performance fetal heart beat detector with color LCD displaying FHR curve, which can meet routine examination requirements of obstetricians. It is applicable for the private obstetrician's office, antepartum clinic, obstetrical department of hospital and remote monitoring situations, or using at home.

### Pictures



### Features

- \* 1.77Inch TFT Color LCD screen
- \* 2Mhz high sensitive probe
- \* Big font FHR display and FHR graphic display
- \* 500MAH rechargeable battery with direct charging function
- \* Real Time FHR and (3 minutes) FHR curve modes

### Specification

#### FHR:

Ultrasound Frequency: 2.0 or 3.0MHz

Intensity : 5mW / cm<sup>2</sup> or less

FHR Range: 60-210bpm

FHR Range for manual count mode: 30-250bpm

FHR Accuracy:  $\pm 1$ bpm