



Fingertip Pulse Oximeter

FDA 510K CE 0123

RoHS



Model: A310



Feature

- Support 6 display modes.
- Dual color OLED display.
- Lightweight and Easy-to-use.
- Visual & sound reminder function.
- Algorithm inside.
- Low perfusion: 0.3%
- Anti – Movement.
- Adjust the parameters in friendly menu.
- Ip22 water splash proof.
- Low battery voltage indicator.
- Automatically switch off while no signal.
- Standard two AAA 1.5V alkaline batteries support more than 20 hours continuous work.
- Optional silicone cover and carry bag.
- With perfusion index display.

Specifications

SpO₂

Measurement Range: 35-100%
 Resolution: 1%
 Measurement Accuracy: $\pm 2\%$ (70%-100%),
 unspecified ($\leq 70\%$)
 Alarm Range: 50% ~ 100%

Pulse Rate:

Measurement Range: 30-250 BPM
 Resolution: 1 BPM
 Measurement Accuracy: ± 2 BPM
 Alarm Range: 30 ~ 250 BPM

Dimension (L×W×H):

60mm×36mm×28mm

Net Weight: 50 g

Color box Dimension (L×W×H):

80mm×60mm×55mm

Gross Weight: 70 g

Outer carton Dimension (L×W×H):

440mm×340mm×290mm

Gross Weight (100pcs): 8.7 Kg

Warranty: 1 Year