



Excellence in Human Service  
and Technology

# Echo Sounder ES-101EX 2MHz FHR Model

## HAND-HELD FETAL DOPPLER WITH LCD

Detects fetal heartbeat and displays the rate.

## ACOUSTIC FEEDBACK CANCELLER

For clear sounds with built-in speaker & transducer.

## EASY & SIMPLE OPERATION

One-handed operation includes freezing heart rate for notation.

## CLINICAL APPLICATIONS

- Detection of fetal heart sounds
- Determination and monitoring of fetal heart rate
- Verification of fetal heart viability by patient trauma

## DISPLAY SAMPLES

Low Battery  
Heart Beat  
Sound



Heart Rate

## SPECIFICATIONS

Power:	DC 3V, 2 AA alkaline batteries
Battery life:	5 hours or more
Automatic shut-off:	3 min. with no signal
Speaker output:	Max. 250 mW
Frequency:	2.25 MHz
Heart rate range:	60 - 240 bpm
LCD display:	Heart rate, sound volume, low battery indicator
Output:	Headset output
Dimensions:	45 (W) x 135 (D) x 30 (H) mm
Weight:	Approx. 150 grams, including batteries
Electrical safety:	Conforms to IEC-60601-1

\*Specifications subject to change.



**Hadeco, Inc.**  
Kawasaki, Japan

www.hadeco.co.jp/

Sole Distributor:  
**Smart Nine, Corp.**  
Kawasaki, Japan

TEL +81-44-860-3901  
FAX +81-44-856-0179



101EX2 CL E S-10  
080-00126-00  
Printed in Japan